

## 兔抗 APOL4 多克隆抗体

中文名称：兔抗 APOL4 多克隆抗体

英文名称： Anti-APOL4 rabbit polyclonal antibody

别名： APOLIV; APOL-IV

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： APOL4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the apolipoprotein L family and may play a role in lipid exchange and transport throughout the body, as well as in reverse cholesterol transport from peripheral cells to the liver. Two transcript variants encoding two different isoforms have been found for this gene. Only one of the isoforms appears to be a secreted protein.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	APOL4
<b>Immunogen:</b>	Synthetic peptide of human APOL4
<b>Full name:</b>	apolipoprotein L4

<b>Synonyms:</b>	APOLIV; APOL-IV
<b>SwissProt:</b>	Q9BPW4
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	10-50
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	HL60 cell lysate
<b>WB Recommended dilution:</b>	200-1000



