

兔抗 APOL6 多克隆抗体

中文名称：兔抗 APOL6 多克隆抗体

英文名称：Anti-APOL6 rabbit polyclonal antibody

别名：APOLVI; APOL-VI

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：APOL6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. Widely expressed; highly expressed in the uterus, fetal brain and spinal cord, also detected in heart, liver, lung, colon, spleen, thymus, prostate, placenta, adrenal gland, salivary and mammary gland.
Applications:	ELISA, IHC
Name of antibody:	APOL6
Immunogen:	Synthetic peptide of human APOL6

Full name:	apolipoprotein L, 6
Synonyms :	APOLVI; APOL-VI
SwissProt:	Q9BWW8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	50-200

