

## 兔抗 APOM 多克隆抗体

中文名称：兔抗 APOM 多克隆抗体

英文名称：Anti-APOM rabbit polyclonal antibody

别名：G3a; NG20; apo-M; HSPC336

相关类别：一抗

抗原：APOM

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is an apolipoprotein and member of the lipocalin protein family. It is found associated with high density lipoproteins and to a lesser extent with low density lipoproteins and triglyceride-rich lipoproteins. The encoded protein is secreted through the plasma membrane but remains membrane-bound, where it is involved in lipid transport. Alternate splicing results in both coding and non-coding variants of this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	APOM

<b>Immunogen:</b>	Fusion protein of human APOM
<b>Full name:</b>	apolipoprotein M
<b>Synonyms :</b>	G3a; NG20; apo-M; HSPC336
<b>SwissProt:</b>	O95445
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	21 kDa
<b>WB Positive control:</b>	Human fetal liver and liver cancer tissue, human placenta tissue
<b>WB Recommended dilution:</b>	200-1000

