

兔抗 CLSTN1 多克隆抗体

中文名称：兔抗 CLSTN1 多克隆抗体

英文名称： Anti-CLSTN1 rabbit polyclonal antibody

别名： calsyntenin 1; CST-1; CSTN1; CDHR12; PIK3CD; ALC-ALPHA; XB31alpha; alalpha1; alalpha2

相关类别： 一抗

抗原： CLSTN1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene is a member of the calsyntenin family, a subset of the cadherin superfamily. The encoded transmembrane protein, also known as alcadein-alpha, is thought to bind to kinesin-1 motors to mediate the axonal anterograde transport of certain types of vesicle. A myloid precursor protein (APP) is trafficked via these vesicles and so this protein is being investigated to see how it might contribute to the mechanisms underlying Alzheimer's disease. Alternative splicing results in multiple transcript variants.

Applications:	ELISA, IHC
Name of antibody:	CLSTN1
Immunogen:	Synthetic peptide of human CLSTN1
Full name:	calsyntenin 1
Synonyms:	CST-1; CSTN1; CDHR12; PIK3CD; ALC-ALPHA; XB31alpha; alcalpha1; alcalpha2
SwissProt:	O94985
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human ovarian cancer
IHC Recommend dilution:	10-50



