

兔抗 AQP4 多克隆抗体

中文名称：兔抗 AQP4 多克隆抗体

英文名称：Anti-AQP4 rabbit polyclonal antibody

别名：MIWC; HMIWC2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AQP4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the aquaporin family of intrinsic membrane proteins that function as water-selective channels in the plasma membranes of many cells. The encoded protein is the predominant aquaporin found in brain. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	AQP4
Immunogen:	Synthetic peptide of human AQP4
Full name:	aquaporin 4

Synonyms :	MIWC; HMIWC2
SwissProt:	P55087
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	35 kDa
WB Positive control:	Mouse kidney tissue
WB Recommended dilution:	200-1000



