

兔抗 AQP11 多克隆抗体

中文名称: 兔抗 AQP11 多克隆抗体

英文名称: Anti-AQP11 rabbit polyclonal antibody

别名: AQPX1

抗原: AQP11

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Aquaporins facilitate the transport of water and small neutral solutes across cell membranes.
Applications:	ELISA, IHC
Name of antibody:	AQP11
Immunogen:	Synthetic peptide of human AQP11
Full name:	aquaporin 11
Synonyms :	AQPX1
SwissProt:	Q8NBQ7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

