

兔抗 CLUL1 多克隆抗体

中文名称：兔抗 CLUL1 多克隆抗体

英文名称：Anti-CLUL1 rabbit polyclonal antibody

别名：clusterin like 1 RA337M

抗原：CLUL1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	CLUL1 is a retina-specific clusterin-like protein which shows nearly 25% identity to clusterin at the protein level. It is expressed specifically in cone photoreceptor cells.
Applications:	ELISA, IHC
Name of antibody:	CLUL1
Immunogen:	Fusion protein of human CLUL1
Full name:	clusterin like 1
Synonyms:	RA337M
SwissProt:	Q15846
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer

IHC Recommend dilution:	20-100
-------------------------	--------

