

## 兔抗 CMC1 多克隆抗体

中文名称：兔抗 CMC1 多克隆抗体

英文名称：Anti-CMC1 rabbit polyclonal antibody

别名：C-X9-C motif containing 1; C3orf68

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CMC1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of the MITRAC (mitochondrial translation regulation assembly intermediate of cytochrome c oxidase complex) complex, that regulates cytochrome c oxidase assembly.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CMC1
<b>Immunogen:</b>	Fusion protein of human CMC1
<b>Full name:</b>	C-X9-C motif containing 1
<b>Synonyms:</b>	C3orf68
<b>SwissProt:</b>	Q7Z7K0
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human esophagus cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	150-300



