

## 兔抗 CMKLR1 多克隆抗体

中文名称: 兔抗 CMKLR1 多克隆抗体

英文名称: Anti-CMKLR1 rabbit polyclonal antibody

别名: DEZ; RVER1; ChemR23; CHEMERINR

抗原: CMKLR1

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

#### Background:

Receptor for the chemoattractant adipokine chemerin/RARRES2 and for the omega-3 fatty acid derived molecule resolvin E1. Interaction with RARRES2 induces activation of intracellular signaling molecules, such as SKY, MAPK1/3 (ERK1/2), MAPK14/P38MAPK and PI3K leading to multifunctional effects, like, reduction of immune responses, enhancing of adipogenesis and angiogenesis. Resolvin E1 down-regulates cytokine production in macrophages by reducing the activation of MAPK1/3 (ERK1/2) and NF-kappa-B. Positively regulates adipogenesis and adipocyte metabolism.

	ism. Acts as a coreceptor for several SIV strains (SIVMAC 316, SIVMAC239, SIVMACL7E-FR and SIVSM62A), as well as a primary HIV-1 strain (92UG024-2).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CMKLR1
<b>Immunogen:</b>	Synthetic peptide of human CMKLR1
<b>Full name:</b>	chemerin chemokine-like receptor 1
<b>Synonyms :</b>	DEZ; RVER1; ChemR23; CHEMERINR
<b>SwissProt:</b>	Q99788
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

