

兔抗 CMTM1 多克隆抗体

中文名称：兔抗 CMTM1 多克隆抗体

英文名称： Anti-CMTM1 rabbit polyclonal antibody

别名： CKLF like MARVEL transmembrane domain containing 1; CKLFH; CKLFH1; CKLFSF1

相关类别： 一抗

抗原： CMTM1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development. Alternatively spliced transcript variants encoding different isoforms have been identified. Naturally occurring read-through transcription occurs between this locus and the neighboring locus CKLF (chemokine-like factor).
Applications:	ELISA, IHC
Name of antibody:	CMTM1

Immunogen:	Synthetic peptide of human CMTM1
Full name:	CKLF like MARVEL transmembrane domain containing 1
Synonyms:	CKLFH; CKLFH1; CKLFSF1
SwissProt:	Q8IZ96
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100

