

兔抗 CMTM2 多克隆抗体

中文名称：兔抗 CMTM2 多克隆抗体

英文名称： Anti-CMTM2 rabbit polyclonal antibody

别名： CKLFSF2

抗原： CMTM2

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the chemokine-like factor gene superfamily, a novel family that links the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development.
Applications:	ELISA, IHC
Name of antibody:	CMTM2
Immunogen:	Fusion protein of human CMTM2
Full name:	CKLF-like MARVEL transmembrane domain containing 2
Synonyms :	CKLFSF2
SwissProt:	Q8TAZ6

ELISA Recommended dilution:	2000-10000
IHC positive control:	Human colon cancer and human esophagus cancer
IHC Recommend dilution:	100-300

