

## 兔抗 ARHGAP1 多克隆抗体

- 中文名称：兔抗 ARHGAP1 多克隆抗体
- 英文名称：Anti-ARHGAP1 rabbit polyclonal antibody
- 别名：RHOGAP; RHOGAP1; CDC42GAP; p50rhoGAP
- 相关类别：一抗
- 抗原：ARHGAP1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a large family of proteins that activate Rho-type guanosine triphosphate (GTP) metabolizing enzymes. The encoded protein contains a SRC homology 3 domain and interacts with Bcl-2-associated protein family members. GTPase activator for the Rho, Rac and Cdc42 proteins, converting them to the putatively inactive GDP-bound state. Cdc42 seems to be the preferred substrate.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ARHGAP1

<b>Immunogen:</b>	Fusion protein of human ARHGAP1
<b>Full name:</b>	Rho GTPase activating protein 1
<b>Synonyms :</b>	RHOGAP; RHOGAP1; CDC42GAP; p50rhoGAP
<b>SwissProt:</b>	Q07960
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	15-50

