

## 兔抗 ARHGAP4 多克隆抗体

中文名称：兔抗 ARHGAP4 多克隆抗体

英文名称：Anti-ARHGAP4 rabbit polyclonal antibody

别名：C1; RGC1; p115; SrGAP4; RhoGAP4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ARHGAP4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a member of the rhoGAP family of proteins which play a role in the regulation of small GTP-binding proteins belonging to the RAS superfamily. The protein encoded by the orthologous gene in rat is localized to the Golgi complex and can redistribute to microtubules. The rat protein stimulates the activity of some Rho GTPases in vitro. Genomic deletions of this gene and a neighboring gene have been found in patients with nephrogenic diabetes insipidus. Multiple transcript variants encoding different isoforms have been found for this gene.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ARHGAP4
<b>Immunogen:</b>	Synthetic peptide of human ARHGAP4
<b>Full name:</b>	Rho GTPase activating protein 4
<b>Synonyms :</b>	C1; RGC1; p115; SrGAP4; RhoGAP4
<b>SwissProt:</b>	P98171
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

