

兔抗 ARHGAP5 多克隆抗体

- 中文名称：兔抗 ARHGAP5 多克隆抗体
- 英文名称：Anti-ARHGAP5 rabbit polyclonal antibody
- 别名：GF12; p190-B; RhoGAP5; p190BRhoGAP
- 相关类别：一抗
- 抗原：ARHGAP5
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. GTPase-activating protein for Rho family members. May play a role in the reduction of the p21rasGTPase-activating potential of RASA1/p120 GAP.
Applications:	ELISA, IHC
Name of antibody:	ARHGAP5
Immunogen:	Fusion protein of human ARHGAP5

Full name:	Rho GTPase activating protein 5
Synonyms :	GFI2; p190-B; RhoGAP5; p190BRhoGAP
SwissProt:	Q13017
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100

