

兔抗 ARHGAP9 多克隆抗体

中文名称：兔抗 ARHGAP9 多克隆抗体

英文名称： Anti-ARHGAP9 rabbit polyclonal antibody

别名： 10C; RGL1

抗原： ARHGAP9

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the Rho-GAP family of GTPase activating proteins. The protein has substantial GAP activity towards several Rho-family GTPases in vitro, converting them to an inactive GDP-bound state. It is implicated in regulating adhesion of hematopoietic cells to the extracellular matrix. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	ARHGAP9
Immunogen:	Synthesized peptide derived from internal of human RHG9.

Full name:	Rho GTPase activating protein 9
Synonyms :	10C; RGL1
SwissProt:	Q9BRR9
WB Predicted band size:	83 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000

