

兔抗 ARHGAP15 多克隆抗体

中文名称：兔抗 ARHGAP15 多克隆抗体

英文名称：Anti-ARHGAP15 rabbit polyclonal antibody

别名：BM046

相关类别：一抗

抗原：ARHGAP15

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	RHO GTPases regulate diverse biologic processes, and their activity is regulated by RHO GTPase-activating proteins (GAPs), such as ARHGAP15. GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Has activity toward RAC1. Overexpression results in an increase in actin stress fibers and cell contraction.
Applications:	ELISA, WB, IHC
Name of antibody:	ARHGAP15
Immunogen:	Fusion protein of human ARHGAP15
Full name:	Rho GTPase activating protein 15

Synonyms :	BM046
SwissProt:	Q53QZ3
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	55 kDa
WB Positive control:	Jurkat and 293T cells
WB Recommended dilution:	200-1000



