

ARL11 抗原（重组蛋白）

中文名称：ARL11 抗原（重组蛋白）

英文名 ARL11 Antigen (Recombinant Protein)

别名 ARLTS1

储存 冷冻（-20℃）

相关类别 抗原

概述

Full length fusion protein

技术规格

Full name:	ADP-ribosylation factor-like 11
Synonyms:	ARLTS1
Swissprot:	Q969Q4
Gene Accession:	BC013150
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers.