

ARL2 抗原(重组蛋白)

中文名称: ARL2 抗原(重组蛋白)

英文名称: ARL2 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Full length fusion protein

技术规格

Full name:	ADP-ribosylation factor-like 2
Synonyms:	ARFL2
Swissprot:	P36404
Gene Accession:	BC002530
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a small GTP-binding protein of the RAS superfamily which functions as an ADP-ribosylation factor (ARF). The encoded protein is one of a functionally distinct group of ARF-like genes. Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP). GTP-binding protein that does not act as an allosteric activator of the cholera toxin catalytic subunit. Regulates formation of new microtubules and centrosome integrity. Prevents the TBCD-induced microtubule destruction.

