

ARL8A 抗原（重组蛋白）

中文名称： ARL8A 抗原（重组蛋白）

英文名称： ARL8A Antigen (Recombinant Protein)

别名： GIE2; ARL10B

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	ADP ribosylation factor like GTPase 8A
Synonyms:	GIE2; ARL10B
Swissprot:	Q96BM9
Gene Accession:	BC015408
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May play a role in lysosomes motility. Alternatively, may play a role in chromosome segregation.