

兔抗 ARL11 多克隆抗体

中文名称：兔抗 ARL11 多克隆抗体

英文名称： Anti-ARL11 rabbit polyclonal antibody

别名： ARLTS1

相关类别： 一抗

抗原： ARL11

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

| | |
|--------------------------|---|
| Background: | This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers. |
| Applications: | ELISA, IHC |
| Name of antibody: | ARL11 |
| Immunogen: | Fusion protein of human ARL11 |
| Full name: | ADP-ribosylation factor-like 11 |
| Synonyms : | ARLTS1 |
| SwissProt: | Q969Q4 |

| | |
|------------------------------------|----------------------|
| ELISA Recommended dilution: | 2000-5000 |
| IHC positive control: | Human thyroid cancer |
| IHC Recommend dilution: | 25-100 |

