

## 兔抗 COG1 多克隆抗体

中文名称：兔抗 COG1 多克隆抗体

英文名称：Anti-COG1 rabbit polyclonal antibody

别名：LDLB; CDG2G

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：COG1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of eight proteins (Cog1-8) which form a Golgi-localized complex (COG) required for normal Golgi morphology and function. It is thought that this protein is required for steps in the normal medial and trans Golgi-associated processing of glycoconjugates and plays a role in the organization of the Golgi-localized complex.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	COG1
<b>Immunogen:</b>	Synthetic peptide of human COG1
<b>Full name:</b>	component of oligomeric golgi complex 1

<b>Synonyms :</b>	LDLB; CDG2G
<b>SwissProt:</b>	Q8WTW3
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

