

## C11orf49 抗原（重组蛋白）

中文名称：C11orf49 抗原（重组蛋白）

英文名称：C11orf49 Antigen (Recombinant Protein)

别名：chromosome 11 open reading frame 49

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-274 amino acids of human C11orf49

技术规格

<b>Full name:</b>	chromosome 11 open reading frame 49
<b>Swissprot:</b>	Q9H6J7
<b>Gene Accession:</b>	BC009455
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The function of this protein remains unknown.