

C4orf46 抗原（重组蛋白）

中文名称：C4orf46 抗原（重组蛋白）

英文名称：C4orf46 Antigen (Recombinant Protein)

别名：RCDG1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	chromosome 4 open reading frame 46
Synonyms:	RCDG1
Swissprot:	Q504U0
Gene Accession:	BC065017
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a small, conserved protein of unknown function that is expressed in a variety of tissues. There are pseudogenes for this gene on chromosomes 6, 8, 16, and X. Alternative splicing results in multiple transcript variants.