

兔抗 COG4 多克隆抗体

- 中文名称：兔抗 COG4 多克隆抗体
- 英文名称：Anti-COG4 rabbit polyclonal antibody
- 别名：component of oligomeric golgi complex 4 COD1; CDG2J
- 相关类别：一抗
- 抗原：COG4
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	COG4
Immunogen:	Fusion protein of human COG4
Full name:	component of oligomeric golgi complex 4
Synonyms:	COD1; CDG2J; SWILS

SwissProt:	Q9H9E3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	100-200



