

## 兔抗 ARPC5L 多克隆抗体

中文名称：兔抗 ARPC5L 多克隆抗体

英文名称： Anti-ARPC5L rabbit polyclonal antibody

别名： actin related protein 2/3 complex subunit 5 like; ARC16-2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ARPC5L

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ARPC5L
<b>Immunogen:</b>	Fusion protein of human ARPC5L
<b>Full name:</b>	actin related protein 2/3 complex subunit 5 like
<b>Synonyms:</b>	ARC16-2
<b>SwissProt:</b>	Q9BPX5

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	30-150

