

兔抗 COL3A1 多克隆抗体

中文名称：兔抗 COL3A1 多克隆抗体

英文名称： Anti-COL3A1 rabbit polyclonal antibody

别名： EDS4A

抗原： COL3A1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the pro-alpha1 chains of type III collagen, a fibrillar collagen that is found in extensible connective tissues such as skin, lung, uterus, intestine and the vascular system, frequently in association with type I collagen. Mutations in this gene are associated with Ehlers-Danlos syndrome types IV, and with aortic and arterial aneurysms. Two transcripts, resulting from the use of alternate polyadenylation signals, have been identified for this gene.
Applications:	ELISA, IHC
Name of antibody:	COL3A1

Immunogen:	Fusion protein of human COL3A1
Full name:	collagen, type III, alpha 1
Synonyms :	EDS4A
SwissProt:	P02461
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human lung cancer and human prostate cancer
IHC Recommend dilution:	50-200

