

## CARTPT 抗原（重组蛋白）

中文名称：CARTPT 抗原（重组蛋白）

英文名称：CARTPT Antigen (Recombinant Protein)

别名：CART

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	CART prepropeptide
<b>Synonyms:</b>	CART
<b>Swissprot:</b>	Q16568
<b>Gene Accession:</b>	BC029882
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a secreted protein which is processed by prohormone/proprotein convertases to produce smaller, biologically active peptides. Expression of the transcript for this gene is regulated by certain drugs such as cocaine, and the encoded protein is thought to be involved in the regulation of appetite and stress. Mutations in this gene are associated with susceptibility to obesity.