

## CDC37L1 抗原（重组蛋白）

中文名称：CDC37L1 抗原（重组蛋白）

英文名称： CDC37L1 Antigen (Recombinant Protein)

别名： HARC; CDC37B; RP11-6J24.5

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	cell division cycle 37-like 1
<b>Synonyms:</b>	HARC; CDC37B; RP11-6J24.5
<b>Swissprot:</b>	Q7L3B6
<b>Gene Accession:</b>	BC014133
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding