

兔抗 ARSE 多克隆抗体

中文名称：兔抗 ARSE 多克隆抗体

英文名称： Anti-ARSE rabbit polyclonal antibody

别名： ASE; CDPX; CDPX1; CDPXR

抗原： ARSE

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Arylsulfatase E is a member of the sulfatase family. It is glycosylated postranslationally and localized to the golgi apparatus. Sulfatases are essential for the correct composition of bone and cartilage matrix. X-linked chondrodysplasia punctata, a disease characterized by abnormalities in cartilage and bone development, has been linked to mutations in this gene. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the Y chromosome.
Applications:	WB

Name of antibody:	ARSE
Immunogen:	Synthesized peptide derived from Internal of human ARSE.
Full name:	arylsulfatase E (chondrodysplasia punctata 1)
Synonyms :	ASE; CDPX; CDPX1; CDPXR
SwissProt:	P51690
WB Predicted band size:	66 kDa
WB Positive control:	K562, 293 and Jurkat cells lysates
WB Recommended dilution:	500-3000

