

## 兔抗 COL12A1 多克隆抗体

中文名称：兔抗 COL12A1 多克隆抗体

英文名称： Anti-COL12A1 rabbit polyclonal antibody

别名： COL12A1L; BA209D8.1; DJ234P15.1

抗原： COL12A1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes the alpha chain of type XII collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XII collagen is a homotrimer found in association with type I collagen, an association that is thought to modify the interactions between collagen I fibrils and the surrounding matrix. Alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Applications:</b>	IHC, IF
<b>Name of antibody:</b>	COL12A1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human Collage

	n XII $\alpha$ 1.
<b>Full name:</b>	collagen, type XII, alpha 1
<b>Synonyms :</b>	COL12A1L; BA209D8.1; DJ234P15.1
<b>SwissProt:</b>	Q99715
<b>IHC positive control:</b>	Human skeletal muscle tissue
<b>IHC Recommend dilution:</b>	50-100
<b>IF positive control:</b>	COS-7 cells
<b>IF Recommend dilution:</b>	100-500



