

兔抗 ARSG 多克隆抗体

中文名称：兔抗 ARSG 多克隆抗体

英文名称：Anti-ARSG rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）

抗原：ARSG

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the sulfatase enzyme family. Sulfatases hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules. This protein displays arylsulfatase activity at acidic pH, as is typical of lysosomal sulfatases, and has been shown to localize in the lysosomes. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ARSG
Immunogen:	Fusion protein of human ARSG

Full name:	arylsulfatase G
SwissProt:	Q96EG1
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer and Human gastric cancer
IHC Recommend dilution:	50-200

