

## 兔抗 COL14A1 多克隆抗体

中文名称：兔抗 COL14A1 多克隆抗体

英文名称： Anti-COL14A1 rabbit polyclonal antibody

别名： UND

抗原： COL14A1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis. [
<b>Applications:</b>	IHC, IF
<b>Name of antibody:</b>	COL14A1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human Collagen XIV $\alpha$ 1.
<b>Full name:</b>	collagen, type XIV, alpha 1
<b>Synonyms :</b>	UND

<b>SwissProt:</b>	Q05707
<b>IHC positive control:</b>	Human heart tissue
<b>IHC Recommend dilution:</b>	50-100
<b>IF positive control:</b>	Hela cells
<b>IF Recommend dilution:</b>	100-500



