

兔抗 COL15A1 多克隆抗体

中文名称: 兔抗 COL15A1 多克隆抗体

英文名称: Anti-COL15A1 rabbit polyclonal antibody

抗 原: COL15A1

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the alpha chain of type XV collagen, a member of the FACIT collagen family (fibril-associated collagens with interrupted helices). Type XV collagen has a wide tissue distribution but the strongest expression is localized to basement membrane zones so it may function to adhere basement membranes to underlying connective tissue stroma. The proteolytically produced C-terminal fragment of type XV collagen is restin, a potentially antiangiogenic protein that is closely related to endostatin. Mouse studies have shown that collagen XV deficiency is associated with muscle and microvessel deterioration.
Applications:	ELISA, IHC
Name of antibody:	COL15A1

Immunogen:	Synthetic peptide of human COL15A1
Full name:	collagen, type XV, alpha 1
SwissProt:	P39059
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100

