

兔抗 ARSK 多克隆抗体

中文名称：兔抗 ARSK 多克隆抗体

英文名称：Anti-ARSK rabbit polyclonal antibody

别名：TSULF

相关类别：一抗

储存：冷冻（-20℃）

抗原：ARSK

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Sulfatases (EC 3.1.5.6), such as ARSK, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules.
Applications:	ELISA, IHC
Name of antibody:	ARSK
Immunogen:	Fusion protein of human ARSK
Full name:	arylsulfatase family, member K
Synonyms:	TSULF
SwissProt:	Q6UWY0

IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	50-200

