

兔抗 ART3 多克隆抗体

中文名称：兔抗 ART3 多克隆抗体

英文名称： Anti-ART3 rabbit polyclonal antibody

别名： ARTC3

相关类别： 一抗

抗原： ART3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an arginine-specific ADP-ribosyltransferase. The encoded protein catalyzes a reversible reaction which modifies proteins by the addition or removal of ADP-ribose to an arginine residue to regulate the function of the modified protein. An ADP-ribosyltransferase pseudogene is located on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ART3
Immunogen:	Fusion protein of human ART3

Full name:	ADP-ribosyltransferase 3
Synonyms :	ARTC3
SwissProt:	Q13508
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

