

兔抗 COL19A1 多克隆抗体

中文名称：兔抗 COL19A1 多克隆抗体

英文名称： Anti-COL19A1 rabbit polyclonal antibody

别名： COL9A1L; D6S228E

抗原： COL19A1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes the alpha chain of type XIX collagen, a member of the FACIT collagen family (fibril-associated collagens with interrupted helices). Although the function of this collagen is not known, other members of this collagen family are found in association with fibril-forming collagens such as type I and II, and serve to maintain the integrity of the extracellular matrix. The transcript produced from this gene has an unusually large 3' UTR which has not been completely sequenced.
Applications:	IHC, IF
Name of antibody:	COL19A1

Immunogen:	Synthesized peptide derived from internal of human Collagen XIX α 1.
Full name:	collagen, type XIX, alpha 1
Synonyms :	COL9A1L; D6S228E
SwissProt:	Q14993
IHC positive control:	Human thymus glan tissue
IHC Recommend dilution:	50-100
IF positive control:	A549 cells
IF Recommend dilution:	100-500



