

兔抗 ART5 多克隆抗体

中文名称：兔抗 ART5 多克隆抗体

英文名称： Anti-ART5 rabbit polyclonal antibody

别名： ARTC5

相关类别： 一抗

抗原： ART5

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the ARG-specific ADP-ribosyltransferase family. Proteins in this family regulate the function of target proteins by attaching ADP-ribose to specific amino acid residues in their target proteins. The mouse homolog lacks a glycosylphosphatidylinositol-anchor signal sequence and is predicted to be a secretory enzyme. Transcript variants with different 5' UTRs, but encoding the same protein have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ART5

Immunogen:	Fusion protein of human ART5
Full name:	ADP-ribosyltransferase 5
Synonyms :	ARTC5
SwissProt:	Q96L15
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	100-300

