

兔抗 COL23A1 多克隆抗体

中文名称：兔抗 COL23A1 多克隆抗体

英文名称： Anti-COL23A1 rabbit polyclonal antibody

抗 原： COL23A1

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	COL23A1 is a member of the transmembrane collagens, a subfamily of the nonfibrillar collagens that contain a single pass hydrophobic transmembrane domain (Banyard et al., 2003 [PubMed 12644459]).[supplied by OMIM, Mar 2008]
Applications:	IF
Name of antibody:	COL23A1
Immunogen:	Synthesized peptide derived from internal of human Collagen XXIII α 1.
Full name:	collagen, type XXIII, alpha 1
Synonyms :	0
SwissProt:	Q86Y22
IF positive control:	HepG2 cells
IF Recommend dilution:	100-500

