

兔抗 ARX 多克隆抗体

中文名称：兔抗 ARX 多克隆抗体

英文名称：Anti-ARX rabbit polyclonal antibody

别名：ISSX; PRYS; CT121; EIEE1; MRX29; MRX32; MRX33; MRX36; MRX38; MRX43; MRX54; MRX76; MRX87; MRXS1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ARX

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:

This gene is a homeobox-containing gene expressed during development. The expressed protein contains two conserved domains, a C-peptide (or aristaless domain) and the prd-like class homeobox domain. It is a member of the group-II aristaless-related protein family whose members are expressed primarily in the central and/or peripheral nervous system. This gene is thought to be involved in CNS development. Mutations in this gene cause X-linked mental retardation and epilepsy.

Applications:	WB
Name of antibody:	ARX
Immunogen:	Synthesized peptide derived from internal of human A RX.
Full name:	aristaless related homeobox
Synonyms :	ISSX; PRTS; CT121; EIEE1; MRX29; MRX32; MRX33; MR X36; MRX38; MRX43; MRX54; MRX76; MRX87; MRXS1
SwissProt:	Q96QS3
WB Predicted band size:	58 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

