

CER1 抗原（重组蛋白）

中文名称：CER1 抗原（重组蛋白）

英文名称：CER1 Antigen (Recombinant Protein)

别名：cerberus 1, DAN family BMP antagonist; DAND4

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 18-267 amino acids of human CER1

技术规格

Full name:	cerberus 1, DAN family BMP antagonist
Synonyms:	DAND4
Swissprot:	O95813
Gene Accession:	BC069491
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a cytokine member of the cysteine knot superfamily, characterized by nine conserved cysteines and a cysteine knot region. The cerberus-related cytokines, together with Dan and DRM/Grem1, represent a group of bone morphogenetic protein (BMP) antagonists that can bind directly to BMPs and inhibit their activity.

mlbio 梅联生物
Good elisakit producers