

兔抗 AS3MT 多克隆抗体

中文名称：兔抗 AS3MT 多克隆抗体

英文名称：Anti-AS3MT rabbit polyclonal antibody

别名：arsenite methyltransferase CYT19

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AS3MT

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	AS3MT catalyzes the transfer of a methyl group from S-adenosyl-L-methionine (AdoMet) to trivalent arsenical and may play a role in arsenic metabolism (Lin et al., 2002 [PubMed 11790780]).
Applications:	ELISA, IHC
Name of antibody:	AS3MT
Immunogen:	Synthetic peptide of human AS3MT
Full name:	arsenite methyltransferase
Synonyms:	CYT19
SwissProt:	Q9HBK9
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100



