

兔抗 COMMD10 多克隆抗体

中文名称：兔抗 COMMD10 多克隆抗体

英文名称：Anti-COMMD10 rabbit polyclonal antibody

别名：PTD002

抗原：COMMD10

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

相关类别：一抗

技术规格

Background:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966).
Applications:	ELISA, WB, IHC
Name of antibody:	COMMD10
Immunogen:	Fusion protein of human COMMD10
Full name:	COMM domain containing 10
Synonyms:	PTD002
SwissProt:	Q9Y6G5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human colorectal cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	23 kDa
WB Positive control:	HEPG2,293T and Hela cell lysates
WB Recommended dilution:	200-1000



