

## CHMP1A 抗原（重组蛋白）

中文名称：CHMP1A 抗原（重组蛋白）

英文名称：CHMP1A Antigen (Recombinant Protein)

别名：PCH8; CHMP1; PRSM1; PCOLN3; VPS46A; VPS46-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	charged multivesicular body protein 1A
<b>Synonyms:</b>	PCH8; CHMP1; PRSM1; PCOLN3; VPS46A; VPS46-1
<b>Swissprot:</b>	Q9HD42
<b>Gene Accession:</b>	BC132711
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the CHMP/Chmp family of proteins which are involved in multivesicular body sorting of proteins to the interiors of lysosomes. The initial prediction of the protein sequence encoded by this gene suggested that the encoded protein was a metalloproteinase. The nomenclature has been updated recently to reflect the correct biological function of this encoded protein. Several transcripts encoding different isoforms have been found for this gene.

e.