

兔抗 ASB2 多克隆抗体

中文名称：兔抗 ASB2 多克隆抗体

英文名称： Anti-ASB2 rabbit polyclonal antibody

别名： ASB-2

相关类别： 一抗

抗原： ASB2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the ankyrin repeat and SOCS box-containing (ASB) protein family. These proteins play a role in protein degradation by coupling suppressor of cytokine signalling (SOCS) proteins with the elongin BC complex. The encoded protein is a subunit of a multimeric E3 ubiquitin ligase complex that mediates the degradation of actin-binding proteins. This gene plays a role in retinoic acid-induced growth inhibition and differentiation of myeloid leukemia cells. Alternatively spliced transcript variants encoding multiple isoforms have been

	observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	ASB2
Immunogen:	Fusion protein of human ASB2
Full name:	ankyrin repeat and SOCS box containing 2
Synonyms :	ASB-2
SwissProt:	Q96Q27
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

