

兔抗 ASB4 多克隆抗体

中文名称：兔抗 ASB4 多克隆抗体

英文名称： Anti-ASB4 rabbit polyclonal antibody

别名： ASB-4

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： ASB4

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene but some of the full length sequences are not known.
Applications:	ELISA, IHC

Name of antibody:	ASB4
Immunogen:	Fusion protein of human ASB4
Full name:	ankyrin repeat and SOCS box containing 4
Synonyms :	ASB-4
SwissProt:	Q9Y574
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

