

## CHRAC1 抗原（重组蛋白）

中文名称： CHRAC1 抗原（重组蛋白）

英文名称： CHRAC1 Antigen (Recombinant Protein)

别名： chromatin accessibility complex 1; YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15; CHRAC-15

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 2-131 amino acids of human CHRAC1

技术规格

<b>Full name:</b>	chromatin accessibility complex 1
<b>Synonyms:</b>	YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15; CHRAC-15
<b>Swissprot:</b>	Q9NRG0
<b>Gene Accession:</b>	BC015891
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	CHRAC1 is a histone-fold protein that interacts with other histone-fold proteins to bind DNA in a sequence-independent manner. These histone-fold protein dimers combine within larger enzymatic complexes for DNA transcription, replication, and packaging.