

## CIDEB 抗原（重组蛋白）

中文名称：CIDEB 抗原（重组蛋白）

英文名称：CIDEB Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	cell death-inducing DFFA-like effector b
<b>Swissprot:</b>	Q9UHD4
<b>Gene Accession:</b>	BC035970
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This proten can activates apoptosis.It is inhibited by DFFB and intera cts with DFFA and DFFB. It has high levels expression in liver and s mall intestine, and has lower levels expression in colon, kidney and s pleen.CideB specifically interacts with VLDL structural protein, apolipo protein B100 (apoB100), but not with albumin, a PTV cargo protein. I t can interacts with COPII components, Sar1 and Sec24