

兔抗 COPS7A 多克隆抗体

中文名称：兔抗 COPS7A 多克隆抗体

英文名称：Anti-COPS7A rabbit polyclonal antibody

别名：COP9 signalosome subunit 7A CSN7; CSN7A; SGN7a

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：COPS7A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described.
Applications:	ELISA, IHC
Name of antibody:	COPS7A
Immunogen:	Fusion protein of human COPS7A
Full name:	COP9 signalosome subunit 7A
Synonyms:	CSN7; CSN7A; SGN7a

SwissProt:	Q9UBW8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100



