

兔抗 ASB10 多克隆抗体

中文名称：兔抗 ASB10 多克隆抗体

英文名称： Anti-ASB10 rabbit polyclonal antibody

别名： ankyrin repeat and SOCS box containing 10; GLC1F

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ASB10

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. The SOCS box serves to couple suppressor of cytokine signaling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	ASB10
Immunogen:	Fusion protein of human ASB10
Full name:	ankyrin repeat and SOCS box containing 10

Synonyms:	GLC1F
SwissProt:	Q8WXI3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	25-50

