

兔抗 CORO1B 多克隆抗体

中文名称：兔抗 CORO1B 多克隆抗体

英文名称：Anti-CORO1B rabbit polyclonal antibody

别名：coronin 1B; CORONIN-2

相关类别：一抗

抗原：CORO1B

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Members of the coronin family, such as CORO1B, are WD repeat-containing actin-binding proteins that regulate cell motility (Cai et al., 2005 [PubMed 16027158]). Regulates leading edge dynamics and cell motility in fibroblasts. May be involved in cytokinesis and signal transduction (By similarity).
Applications:	ELISA, IHC
Name of antibody:	CORO1B
Immunogen:	Synthetic peptide of human CORO1B
Full name:	coronin 1B
Synonyms:	CORONIN-2

SwissProt:	Q9BR76
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100



